

Microscope Design, Operation, and Alignment

Goal

You will be introduced to the technical construction and the main components of a transmission electron microscope (Philips CM20), the basic operation of this instrument, alignment procedures, image recording, and darkroom procedures. The specimen of this laboratory is a carbon film supporting MoO_3 crystallites. The morphology of these crystallites uniquely correlates with the orientation of the crystal lattice. This is useful for calibrating the rotation of the image and the diffraction pattern as a function of magnification.

Experiment

1. Record pictures of the instrument and specimen holder with a digital camera.
2. Insert a specimen and start up the microscope.
3. Familiarize yourself with the basic operation (apertures, condenser deflection, specimen translation, magnification, focus).
4. Perform a basic alignment (eucentric height, gun tilt, gun shift, condenser aperture).
5. Observe the rotation of the image as a function of magnification.
6. Observe diffraction patterns.
7. Choose a fixed camera length for diffraction patterns and employ double exposures to document the rotation of the image versus the diffraction pattern as a function of image magnification.
8. Develop the negatives and reload the microscope camera.

Report

- Referring to your digital pictures, summarize what you learned about the construction of a transmission electron microscope.
- Summarize the purpose and the individual steps of the alignment procedure.
- Include a plot of image rotation versus the diffraction pattern as a function of magnification and comment on the underlying physics.